

Le Mieux Regenerating Facial with EpiPin

1. Cleanse

Cleanse with appropriate cleanser depending upon skin type. Rinse well and pat dry.

2. Exfoliate

Combine Hungarian Enzyme powder with warm water at a 1:1.5 ratio. Use a small fan brush to apply to the face and neck. Gently massage for a minute. Let dry 5-10 min then remove with warm moist towel. **Tip** - you may use ultrasonic spatula (Skin Perfecter tool) to gently exfoliate using wet 2x2's, one pass.

3. EpiPin

EpiPin on dry skin, then apply Hyaluronic Serum to small sections at a time (4x4 area). You may optionally add a few drops of Vita-C Serum (balances uneven skintone and collagen synthesis), or Rx Complex Serum (great for expression lines and wrinkles) and EpiPin again. Massage in any excess serum into the skin.

4. Mask

Apply Bio Cell+ Mask. Remove paper-like material from mask and spread gel like sheet mask on to face. **Option:** cold hammer, ice globes or ice roller may be glided over mask. Remove mask after 15–20 minutes. Do not rinse.

5. Finish

Massage solution from mask gently into skin. Spritz ISO Cell Recovery solution on skin. Apply O2 Calming Gel (to calm irritated skin) or Apply a few drops of EGF-DNA or StemCell EGF to skin. **Home care instructions:** Tell client to not wash the face for 12 hours. They may however splash the skin with cool water and carefully pat dry. The skin may be slightly pink for the first few hours. No steam, hot yoga or any activity that produces sweat for 24 hours. This may irritate skin. For the first 24 hours to use tepid water to wash face and use gentle cleanser. They may spritz liberally with the Iso-cell recovery solution and apply a few drops of the Derma Relief Serum or Hyaluronic serum with EGF of choice first day or two after procedure. Client should avoid direct sun exposure and wear SPF of choice. Days 2 to 4 they may experience dry skin as they may have a very gentle micro – peel which is essentially an accelerated desquamation of the stratum corneum. Days 5–7 they will see tighter, more evenly toned and textured skin and start seeing the instant results of the EpiPin Treatment.

For ultimate results, do this treatment every week for four weeks then every two weeks for the last two treatments for a total of six treatments.

* EpiPin is 'Product Neutral' and does not promote or endorse any specific product line. We haven't necessarily performed the protocol and ask you to rely on your skin care expertise to determine what products are best to use on your clients.